SONY®

4-591-742-07 (1)


2021-12

Wireless Print System 无线打印组件

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EU: Sony Europe B.V.
Da Vincilaan 7-D1, 1930 Zaventem, Belgium
UK: Sony Europe B.V.
The Heights, Brooklands, Weybridge,
Surrey KT13 0XW, United Kingdom
CH: Sony Europe B.V., Hoofddorp,
Schlieren/Switzerland Branch
Wiesenstrasse 5, 8952 Schlieren, Switzerland

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Before operating the unit, please read this manual thoroughly and retain it for future reference.

Indications for Use/Intended Use

The Sony Wireless USB (WUSB) device is a medical grade, paired transmitter and receiver, and an optional accessory intended for delivery of medical images over the ultra-wide band from ultrasound, mobile C-ARM, endoscopy and cardiac catheterization laboratory equipment, and other medical image capturing devices for printing on Sony medical printers, which are not intended for diagnostic use. The WUSB is a non-sterile reusable device intended for use by qualified healthcare personnel within specified distance between the printer and medical imaging device.

Notes

- This equipment is for medical professionals.
- This equipment is intended for use in medical environments, such as clinics, examination rooms, and operating rooms.

WARNING

To reduce the risk of fire or electric shock, do not expose this equipment to rain or moisture.

To avoid electrical shock, do not open the cabinet. Refer servicing to qualified personnel only.

No modification of this equipment is allowed.

Symbols on the product

This symbol indicates the radio emission equipment.

Consult the instructions for use

Follow the directions in the instructions for use for parts of the unit on which this symbol appears.

This symbol indicates the manufacturer, and appears next to the manufacturer's name and address.

This symbol indicates the Importer, and appears next to the Importer's name and registered office address.

This symbol indicates the European Community representative, and appears next to the European Community representative's name and address.

This symbol indicates the UK Responsible Person, and appears next to the UK Responsible Person's name and address.

This symbol indicates the Swiss authorized representative, and appears next to the Swiss authorized representative's name and address.

This symbol indicates the medical device in the European Community.

This symbol indicates the date of manufacture.

This symbol indicates the serial number.

This symbol indicates the Unique Device Identifier (UDI), and appears next to the bar code representation of the Unique Device Identification.

Storage and transport temperature

This symbol indicates the acceptable temperature range for storage and transport environments.

Storage and transport humidity

This symbol indicates the acceptable humidity range for storage and transport environments.

Storage and transport pressure

This symbol indicates the acceptable atmospheric pressure range for storage and transport environments.

Warning

Keep this product at least 22 cm away from cardiac pacemakers. Radio waves from this product may adversely affect the operation of cardiac pacemakers.

If a fault occurs during operation of this product, stop use of this product.

You are cautioned that any changes or modifications not expressly approved in this manual could void your authority to operate this equipment.

Caution

If the unit malfunctions, do not touch any unit circuitry and the patient at the same time. It may generate a voltage that could be

harmful to the patient. Stop using the unit immediately, and disconnect and remove the unit.

Caution

If a malfunction occurs, such as a printing error, heat build-up may occur on the surface of the UPA-WU10T due to UPA-WU10T circuit failure and may result in a burn injury.

Turn off the medical equipment connected to the UPA-WU10T and wait until the UPA-WU10T completely cools down.

Caution

For cautions and further information on using this product, refer to the instructions for use for the printer.

UPA-WU10R only

IMPORTANT

The nameplate is located on the bottom.

For the customers in the U.S.A.

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications to this equipment. Such modifications could void the user's authority to operate the equipment.

UWB systems operating under the provisions of this section shall bear the following or similar statement in a

conspicuous location on the device or in the instruction manual supplied with the device:

“This equipment may only be operated indoors. Operation outdoors is in violation of 47 U.S.C. 301 and could subject the operator to serious legal penalties.”

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines.

This equipment has very low levels of RF energy that is deemed to comply without testing of specific absorption rate (SAR).

Caution

Federal law (United States of America) restricts this device to sale by or on the order of a licensed healthcare practitioner.

Rx
ONLY

For customers in Canada

This equipment complies with IC radiation exposure limits set forth for an uncontrolled environment and meets RSS-102 of the IC radio frequency (RF) Exposure rules.

This equipment has very low levels of RF energy that is deemed to comply without testing of specific absorption rate (SAR).

This device complies with Industry Canada's licence-exempt RSSs. Operation is subject to the following two conditions:

- (1) This device may not cause interference; and
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

The term “IC:” before the radio certification number only signifies that Industry Canada technical specifications were met.

For customers in Europe

Hereby, Sony Corporation declares that the radio equipment type UPA-WU10T and UPA-WU10R are in compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address:

<https://compliance.sony.eu>

Hereby, Sony Corporation declares that the radio equipment type UPA-WU10T and UPA-WU10R are in compliance with the UK relevant statutory requirements.

The full text of the declaration of conformity is available at the following internet address:

<https://compliance.sony.co.uk>

This equipment is restricted to indoor operations only.

These restrictions apply in the following countries:

BE/BG/CZ/DK/DE/EE/IE/EL/ES/FR/HR/IT/CY/LV/LT/LU/HU/MT/NL/AT/PL/PT/RO/SI/SK/FI/SE/UK/UK(NI)/IS/LI/NO/CH.

Important EMC notices for use in medical environments

- The UPA-WU10T/UPA-WU10R needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the instructions for use.
- The UPA-WU10T/UPA-WU10R is intended for use in a professional healthcare facility environment.
- Portable and mobile RF communications equipment, such as cellular phones, can affect the UPA-WU10T/UPA-WU10R.

Warning

- Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the UPA-WU10T/UPA-WU10R. Otherwise, degradation of the performance of this equipment could result.
- If the UPA-WU10T/UPA-WU10R will be used adjacent to or stacked with other equipment, normal operation of the UPA-WU10T/UPA-WU10R under such configurations should be verified via observation.
- The use of accessories and cables other than those specified, with the exception of replacement parts sold by Sony Corporation, may result in increased emissions or decreased immunity of the UPA-WU10T/UPA-WU10R.

List of cables used for EMC test	
Type of cable	Specifications
USB cables (1-757-429-21) (2)	1 m, Shielded cable

Guidance and manufacturer’s declaration – electromagnetic emissions		
The UPA-WU10T/UPA-WU10R is intended for use in the electromagnetic environment specified below. The customer or the user of the UPA-WU10T/UPA-WU10R should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The UPA-WU10T/UPA-WU10R uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The UPA-WU10T/UPA-WU10R is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity			
The UPA-WU10T/UPA-WU10R is intended for use in the electromagnetic environment specified below. The customer or the user of the UPA-WU10T/UPA-WU10R should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD)	±8 kV contact	±8 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, a relative humidity of at least 30% is recommended.
IEC 61000-4-2	±15 kV air	±15 kV air	
Electrical fast transient/burst	±2 kV for power supply lines	±2 kV for power supply lines	
IEC 61000-4-4	±1 kV for input/output lines	±1 kV for input/output lines	
Surge	±1 kV line(s) to line(s)	±1 kV differential mode	
IEC 61000-4-5	±2 kV line(s) to earth	±2 kV common mode	
Voltage dips, short interruptions and voltage variations on power supply input lines	0% U_T (100% dip in U_T) for 0.5/1 cycles ^a	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the UPA-WU10T/UPA-WU10R requires continued operation during power mains interruptions, it is recommended that the UPA-WU10T/UPA-WU10R be powered from an uninterruptible power supply or a battery.
IEC 61000-4-11	40% U_T (60% dip in U_T) for 5 cycles		
	70% U_T (30% dip in U_T) for 25/30 cycles ^a (for 0.5 sec)		
	0% U_T (100% dip in U_T) for 250/300 cycles ^a (for 5 sec)		

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			
a For example, 10/12 means 10 cycles at 50 Hz or 12 cycles at 60 Hz.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The UPA-WU10T/UPA-WU10R is intended for use in the electromagnetic environment specified below. The customer or the user of the UPA-WU10T/UPA-WU10R should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands ^c	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the UPA-WU10T/UPA-WU10R, including cables, than the recommended separation distance calculated from the equation appliance to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$
	6 Vrms 150 kHz to 80 MHz in ISM bands ^c	6 Vrms	

Radiated RF	3 V/m	3 V/m	IEC 60601-1-2: 2007
IEC 61000-4-3	80 MHz to 2.7 GHz		$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz IEC 60601-1-2: 2014 $d = 2.0 \sqrt{P}$ 80 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the UPA-WU10T/UPA-WU10R is used exceeds the applicable RF compliance level above, the UPA-WU10T/UPA-WU10R should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the UPA-WU10T/UPA-WU10R.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.
- c The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

Recommended separation distances between portable and mobile RF communications equipment and the UPA-WU10T/UPA-WU10R					
The UPA-WU10T/UPA-WU10R is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the UPA-WU10T/UPA-WU10R can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the UPA-WU10T/UPA-WU10R as recommended below, according to the maximum output power of the communications equipment.					
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m				
	IEC 60601-1-2 : 2007			IEC 60601-1-2 : 2014	
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 2.7 GHz $d = 2.0 \sqrt{P}$
0.01	0.12	0.12	0.23	0.12	0.20
0.1	0.38	0.38	0.73	0.38	0.63
1	1.2	1.2	2.3	1.2	2.0
10	3.8	3.8	7.3	3.8	6.3
100	12	12	23	12	20
For transmitters rated a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.					
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.					
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.					

Guidance and manufacturer's declaration – electromagnetic immunity					
The UPA-WU10T/UPA-WU10R is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the UPA-WU10T/UPA-WU10R. Otherwise, degradation of the performance of this equipment could result.					
Immunity test	Band ^a	Service ^a	Modulation	IEC 60601 test level	Compliance level
Proximity fields from RF wireless communications equipment IEC 61000-4-3	380 – 390 MHz	TETRA 400	Pulse modulation 18 Hz	27 V/m	27 V/m
	430 – 470 MHz	GMRS 460 FRS 460	FM ±5 kHz deviation 1 kHz sine	28 V/m	28 V/m
	704 – 787 MHz	LTE Band 13, 17	Pulse modulation 217 Hz	9 V/m	9 V/m
	800 – 960 MHz	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28 V/m	28 V/m
	1,700 – 1,990 MHz	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation 217 Hz	28 V/m	28 V/m
	2,400 – 2,570 MHz	Bluetooth WLAN 802. 11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28 V/m	28 V/m
	5,100 – 5,800 MHz	WLAN 802. 11 a/n	Pulse modulation 217 Hz	9 V/m	9 V/m
NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.					
a For some services, only the uplink frequencies are included.					

Caution

When you dispose of the unit or accessories, you must obey the laws in the relative area or country and the regulations in the relative hospital regarding environmental pollution.

Caution

Do not use the device in a MR (Magnetic Resonance) environment.

It may cause a malfunction, fire, and unwanted movement.

For the customers in the U.S.A.

SONY LIMITED WARRANTY - Please visit

<http://www.sony.com/psa/warranty> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

For the customers in Canada

SONY LIMITED WARRANTY - Please visit

<http://www.sonybiz.ca/pro/lang/en/ca/article/resources-warranty> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

For the customers in Europe

Sony Professional Solutions Europe - Standard Warranty and Exceptions on Standard Warranty.

Please visit <https://pro.sony/support-services/primesupport/support-professional-solutions-europe-standard-product-warranty> for important information and complete terms and conditions.

For the customers in Korea

SONY LIMITED WARRANTY - Please visit

<http://bpeng.sony.co.kr/handler/BPAS-Start> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

Usage Precautions

On Condensation

If the unit is suddenly taken from a cold to a warm location, or if ambient temperature suddenly rises, moisture may form on the outer surface of the unit and/or inside of the unit. This is known as condensation. If condensation occurs, turn off the unit and wait until the condensation clears before operating the unit. Operating the unit while condensation is present may damage the unit.

On wireless communications

- The communications speed may be reduced or communications may be interrupted, depending on the surrounding environment and obstacles, device configuration, RF signal strength, line congestion, or other causes. This is not a malfunction.
- The supported frequency band of the UWB wireless function of the unit may be used by other wireless devices. If other wireless devices are affected by radio interference, cease using the unit.

Use with electrosurgical knives and similar devices

If this unit is used together with an electrosurgical knife, etc., the unit may not operate normally as a result of adverse effects caused by strong radio waves or voltages from the device. This is not a malfunction.

When you use this unit simultaneously with a device from which strong radio waves or voltages are emitted, confirm the effect of

this before using such devices, and install this unit in a way that minimizes the effect of radio wave interference.

On Security

SONY WILL NOT BE LIABLE FOR DAMAGES OF ANY KIND RESULTING FROM A FAILURE TO IMPLEMENT PROPER SECURITY MEASURES ON TRANSMISSION DEVICES, UNAVOIDABLE DATA LEAKS RESULTING FROM TRANSMISSION SPECIFICATIONS, OR SECURITY PROBLEMS OF ANY KIND.


Overview

The UPA-WU10 Wireless Print System (hereinafter, “system”) is an accessory for Sony printers, comprising a UPA-WU10T transmitter and UPA-WU10R receiver. The system connects medical-grade equipment to a printer wirelessly, replacing a USB cable, enabling the printer to be installed in a location separate from the medical-grade equipment.

Name and Function of Parts

Wireless Print Receiver (UPA-WU10R)

Front panel

① **SIGNAL indicator**

Indicates the RF signal strength.
Ensure that two or more indicator bars are lit for stable communication.


② **LINK indicator**

Lights up when a communication link is established with the transmitter.

③ **Channel indicator**

The specified communication channel indicator lights up.

Rear panel

① CH button

Switch the communication channel.

② Slide switch

Switch between run mode (RUN) and setup mode (SETUP).

③ PRINTER connector (USB type A)

Connect to a Sony printer using the supplied USB cable.

④ (DC IN) connector (USB type B)


Connect to a Sony printer port that supplies USB bus power using the supplied USB cable.

Caution

Do not connect to a printer other than a Sony printer. Doing so may pose a risk of fire and/or electric shock.

Wireless Print Transmitter (UPA-WU10T)

① LINK indicator

Lights up when a communication link is established with the receiver.

② Channel indicator

The specified communication channel indicator lights up.

③ USB connector (USB Type A)

Connect to medical equipment.

Connecting

Before using the system, install the driver software for the printer to connect in the medical equipment. Then, use the following procedure to connect wirelessly.

- 1** Slide the slide switch on the rear panel of the receiver to the RUN position.
- 2** Connect the printer to the PRINTER connector on the rear panel of the receiver using the supplied USB cable.
- 3** Connect the printer port that supplies USB bus power to the DC IN connector on the rear panel of the receiver using the supplied USB cable.
- 4** Connect the transmitter to the USB port of the medical device.
- 5** Turn the printer on.
The LINK indicator on the front panel of the receiver starts flashing.
- 6** Turn the medical device on.
The LINK indicator on the front panel of the transmitter starts flashing.

The receiver and transmitter are connected automatically. When the LINK indicators on both the receiver and transmitter are lit solid, printing from the medical device is enabled.

Notes

- Communication may become unstable if the system is housed within a metal case, installed in close proximity to metal surfaces, or if the RF range becomes shortened.

- For stable communication, ensure that two or more SIGNAL indicator bars are lit on the front panel of the receiver.


Using Multiple Wireless Print Systems

To use multiple wireless print systems within the same area, set the systems to different channels. The receiver and transmitter set to the same channel will connect automatically. Up to three systems can be used at the same time.

Use the following procedure to set the channel. The same channel is set on the transmitter and receiver.

Note

When setting the channel, check that there are no other wireless print systems operating in the vicinity.

- 1** Slide the slide switch on the rear panel of the receiver to the SETUP position.

2 Connect the transmitter to the PRINTER connector on the rear panel of the receiver.

3 Connect the printer port that supplies USB bus power to the DC IN connector on the rear panel of the receiver using the supplied USB cable.

The LINK indicator (**A**) on the front panel of the transmitter starts flashing and the currently specified channel indicator lights up solid (channel 1 by factory default).

4 Press the CH button on the rear panel of the receiver to change channel.

The channel changes and the internal settings are modified each time you press the CH button. The LINK indicator turns off while the internal settings are modified, and the channel indicator flashes three times. After the settings are modified, the changed channel indicator is lit solid and the LINK indicator starts flashing.

Note

Do not disconnect the USB bus power while the LINK indicator is off. The system may no longer operate if the settings information is corrupted.

5 Disconnect the USB bus power and the transmitter from the receiver, in that order.

6 Slide the slide switch on the rear panel of the receiver to the RUN position.

Follow the procedure in the “Connecting” section to connect the printer and medical device.

Notes

- If the channels of the transmitter and receiver are different, reconfigure the

channel settings using the procedure in the “Using Multiple Wireless Print Systems” section.

- This system automatically connects a receiver and transmitter set to the same channel. If there are multiple transmitters and receivers set to the same channel within the same communications area, an unintended combination of devices may be connected, causing printing by an unintended printer. When using multiple devices simultaneously, always set the devices to different channels.

Cleaning

Wiping off dirt or stains on the receiver or transmitter using benzine, thinner, acidic cleaning solution, alkaline cleaning solution, abrasive cleaning solution, or chemically-treated cloth may damage the surface finish. Observe the following when cleaning.

- Wipe surfaces clean using a solution containing 50 to 70 v/v% concentration isopropyl alcohol or 76.9 to 81.4 v/v% concentration ethanol.
- For stubborn dirt or stains, wipe with a soft cloth, such as a cleaning cloth, lightly moistened with neutral detergent diluted with water, and then wipe clean using the isopropyl alcohol solution or ethanol solution above.
- Do not wipe vigorously if there is any dirt adhering to the cloth. Doing so may scratch the surfaces.
- Do not allow any rubber or plastic products to touch the surfaces for extended periods. Do not alter, paint, or scratch the surface finish.

Troubleshooting

If a problem occurs, check the following items. If the problem persists, contact the point of purchase.

Symptom	Check
The LINK indicator does not turn on at all.	• Check that the device that supplies power is also turned on. • Check that the slide switch on the receiver is in the RUN position.
The LINK indicator flashes continuously.	• Check that the transmitter and receiver are set to the same channel. • Check that there is no other transmitter or receiver set to the same channel in the vicinity. • Check that the transmitter and receiver are not out of range.
Communication fails.	• Check that at least two SIGNAL indicator bars are lit. → Try decreasing the separation distance between transmitter and receiver.
Printing is slow to start.	• Check that the LINK indicator is not flashing rapidly for long intervals. → There may be interference due to reflected radio waves. Change the position and orientation of the receiver.

Symptom	Check
Does not connect to the intended printer.	- Check that there is no other system set to the same channel in the vicinity. → Change the channel settings.

Specifications

Note

As of July 2019, this system can be used in the EU and EFTA, the United States, Canada, Japan, and China.

- Communication system
 - UWB (Ultra Wide Band)
- Frequency band
 - 7,392 MHz to 8,448 MHz
 - Band Group #6 (Band #9, #10)
- Maximum range
 - Approx. 10 m (33 ft.)¹⁾ line-of-sight
- Maximum number of devices in use simultaneously
 - 3
- Interface
 - Hi-Speed USB (USB 2.0 compliant)
- Power requirement
 - 5 V DC / 0.5 A (USB bus power)
- Center Frequency
 - 7,656 MHz and 8,184 MHz
- Band Width
 - 528 MHz
- Maximum mean power spectral density
 - −41.3 dBm/MHz
- Maximum peak power
 - 0 dBm
- Dimensions (W×H×D)
 - Transmitter: Approx. 19.5 × 8.5 × 41.3 mm (²⁵/₃₂ × ¹¹/₃₂ × 1 ¹¹/₁₆ in.)
 - Approx. 19.5 × 8.5 × 53.3 mm (²⁵/₃₂ × ¹¹/₃₂ × 2 ¹/₈ in.)
 - (including USB connector)
 - Receiver: Approx. 27.6 × 66.1 × 66.1 mm (1 ¹/₈ × 2 ⁵/₈ × 2 ⁵/₈ in.)
- Mass
 - Transmitter: Approx. 8 g (0.28 oz.)
 - Receiver: Approx. 57 g (2.0 oz.)

Operating Temperature
5 °C to 40 °C (41 °F to 104 °F)

Operating Humidity
20% to 80% (no condensation
allowed)

Operating Pressure
700 hPa to 1,060 hPa

Storage and Transport Temperature
–20 °C to +60 °C (–4 °F to
+140 °F)

Storage and Transport Humidity
20% to 80% (no condensation
allowed)

Storage and Transport Pressure
700 hPa to 1,060 hPa

Supplied accessories
USB cables (1-757-429-21) (2)
Stand (1)
Instructions for Use (1)
Service Contact List (1)
Information for Customers in
Europe (1)

1) The communication speed will vary depending on the distance and obstacles between communication equipment, equipment configuration, RF signal strength, line congestion, and the equipment used. Communication may also fail depending on the RF signal strength.

Design and specifications are subject to change without notice.

Medical Specifications

Protection against harmful ingress of water:
Ordinary

Degree of safety in the presence of flammable anesthetics mixtures with air or with oxygen or nitrous oxide:
Not suitable for use in the presence of flammable anesthetics mixture with air or with oxygen or nitrous oxide

Mode of operation:
Continuous

Notes

- Always verify that the unit is operating properly before use. SONY WILL NOT BE LIABLE FOR DAMAGES OF ANY KIND INCLUDING, BUT NOT LIMITED TO, COMPENSATION OR REIMBURSEMENT ON ACCOUNT OF THE LOSS OF PRESENT OR PROSPECTIVE PROFITS DUE TO FAILURE OF THIS UNIT, EITHER DURING THE WARRANTY PERIOD OR AFTER EXPIRATION OF THE WARRANTY, OR FOR ANY OTHER REASON WHATSOEVER.
- SONY WILL NOT BE LIABLE FOR CLAIMS OF ANY KIND MADE BY USERS OF THIS UNIT OR MADE BY THIRD PARTIES.
- SONY WILL NOT BE LIABLE FOR THE TERMINATION OR DISCONTINUATION OF ANY SERVICES RELATED TO THIS UNIT THAT MAY RESULT DUE TO CIRCUMSTANCES OF ANY KIND.

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© 2016 Sony Corporation