

Baaske GL medical keyboard



Perfectly suited for all areas in which hygiene is a must.

The e-medical keyboard GL01 has been designed for a long service life in the field of medicine. It is characterised by a highquality design and a well thought-out material combination: The front is made of hardened glass, the rest of the housing is made of ABS plastic. The Baaske keyboard GL01 therefore provides important prerequisites for easy and quick cleaning. This optimal disinfection combined with the touch functionality and a switchable touchpad make it an ideal solution in hospitals or doctor's surgeries.

Benefits:

- High quality design
- Alphanumeric glass keyboard
- IP 67 - front sealing
- Can be operated with gloves
- Easy to clean
- Integrated touchpad
- USB port

Baaske keyboard GL01 for sensitive areas where hygiene is a must.



The Baaske keyboard GL01 consists of hardened glass on the front and a robust ABS housing.



The USB port, via which the keyboard can be operated and charged, is provided with a cover so that the medical Baaske keyboard GL01 is easy to clean both on the smooth glass top and all around.



The alphanumeric keyboard has an integrated „Macbook“ touchpad and can be operated with gloves.

Technical data

Baaske keyboard GL01

Number of keys	105 keys
Degree of protection	IP 67
Compatible OS	Microsoft® Windows, Mac OSX, Android
Dimensions (LxWxH)	436 x 160 x 20 mm
Material	Hardened glass, robust ABS housing
Coating	insensitive for fingerprints
Touchpad size (LxW)	119 x 102 mm
Weight	0,87 kg
Color	white
Layout	DE
Interface	USB, Bluetooth & RF via Dongle
Operating temperature	20°C to 60°C
Standards	FCC, CE, ROHS
Scope of delivery	Micro-USB-Kable, cleaning cloth
Item No.	2013126